


## TC MAYO DISSECTING SCISSORS

straight; blunt/blunt

23.0cm (9.2")




**SKU:** 011-100-23

**Categories:** [Open Surgery](#), [Scissors](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                   |
| <b>TYPE</b>           | <a href="#">Dissecting Scissors</a>                    |
| <b>Aesculap</b>       | BC558  |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                   |
| <b>KLS Martin</b>     | 11-910-23  |
| <b>AUTHORS</b>        | <a href="#">MAYO</a>                                   |
| <b>Medicon</b>        | 03.50.73   |
| <b>SPECIFICATIONS</b> | <a href="#">blunt/blunt</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                     |
|-----------------|-------------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>                |
| <b>TYPE</b>     | <a href="#">Dissecting Scissors</a> |
| <b>Aesculap</b> | BC558                               |

|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>  |
| <b>KLS Martin</b>     | 11-910-23                             |
| <b>AUTHORS</b>        | <a href="#">MAYO</a>                  |
| <b>Medicon</b>        | 03.50.73                              |
| <b>SPECIFICATIONS</b> | <a href="#">blunt/blunt, straight</a> |