

## WORKING ELEMENT;

passive cutting; closed handle

for Ø 4mm telescopes; HF-connection up

for mono- and bipolar use; ATTENTION: our bipolar resectoscopy system can be used only with HF-units which have the possibility of bipolar cut in saline solution!



**SKU:** 630-235-60B

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>