


WORKING ELEMENT TITANIUM; BIPOLAR

passive cutting; closed handle

Quick-Lock System; for 0° telescope

ATTENTION: our bipolar resectoscopy system can be used only with HF-units which have the possibility of bipolar cut in saline solution!!



SKU: 582-150-00GUB

Categories: [Endoscopy](#), [GUBBINI/Gyn](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Ellipse , Original
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Ellipse , Original
TYPE	Working Elements