

## 2- LUMEN RETRIEVAL BASKET, STAINLESS

hexagonal, 20 x 40mm; 4 wires

3.3mm; 210cm; for guide wire 0.035"



**SKU:** PK-1FHPK3320230K-2L

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)