

2- LUMEN RETRIEVAL BASKET, STAINLESS

hexagonal, 20 x 40mm; 4 wires

3.3mm; 210cm; for guide wire 0.035"



SKU: PK-1FHPK3320230K-2L

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)