

3-RING SYRINGE WITH GLASS BARREL

Luer-Lock connection; 5 ml



SKU: 087-055-05

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-05
Medicon	09.78.53

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-05
Medicon	09.78.53