

3-RING SYRINGE WITH GLASS BARREL

Luer-Lock connection; 10 ml

SKU: 087-055-10

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-10
Medicon	09.78.54

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-10
Medicon	09.78.54