

3-RING SYRINGE WITH GLASS BARREL

Luer-Lock connection; 2 ml

SKU: 087-055-02

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-02
Medicon	09.78.52

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Syringes
KLS Martin	18-915-02
Medicon	09.78.52