

3-WAY STOPCOCK

1X LL= male; 2X LL= female



SKU: 087-922-00

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Stopcocks and Connecting Pieces |
| Aesculap | SK320 |
| KLS Martin | 18-845-00 |
| Medicon | 09.79.65 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Stopcocks and Connecting Pieces |
| Aesculap | SK320 |
| KLS Martin | 18-845-00 |
| Medicon | 09.79.65 |