

90° DOUBLENEPHROSCOPE

19Fr., working channel Ø 5.0mm

view 15°, working length 200mm

SKU: 601-130-90

Categories: [Endoscopy](#), [Equipment](#), [Urology](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)