

## ADAPTER PTFE TEFLON Ø 43 MM (OUTER)

white; for use with ultrasound



**SKU:** BLUO/02-00045-26

**Categories:** [Blue endo](#), [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uterusmanipulator](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uterusmanipulator](#)