


ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth

16.5cm (6.6")



SKU: 070-860-16

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| Aesculap | BV995 |
| KLS Martin | 15-771-16 |
| AUTHORS | ADSON |
| Medicon | 24.30.16 |
| SPECIFICATIONS | sharp , toothed |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| Aesculap | BV995 |
| KLS Martin | 15-771-16 |

AUTHORS

[ADSON](#)

Medicon

24.30.16

SPECIFICATIONS

[sharp, toothed](#)