

ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth

20.0cm (8")




SKU: 070-860-20

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-770-20
AUTHORS	ADSON
Medicon	24.30.20
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-770-20
AUTHORS	ADSON
Medicon	24.30.20

SPECIFICATIONS

[sharp, toothed](#)