

## ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

16.5cm (6.6")




**SKU:** 070-861-16

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>