

ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

16.5cm (6.6")




SKU: 070-861-16

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed