

## ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

20.0cm (8")




**SKU:** 070-861-20

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-771-20
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	24.30.21
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-771-20
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	24.30.21

## SPECIFICATIONS

[blunt](#), [toothed](#)