

ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

20.0cm (8")




SKU: 070-861-20

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-20
AUTHORS	ADSON
Medicon	24.30.21
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-20
AUTHORS	ADSON
Medicon	24.30.21

SPECIFICATIONS

[blunt](#), [toothed](#)