

ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

26.0cm (10.4")




SKU: 070-861-26

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed