

ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth

13.0cm (5.2")




SKU: 070-860-13

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV085
KLS Martin	15-770-13
AUTHORS	ADSON
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV085
KLS Martin	15-770-13
AUTHORS	ADSON

SPECIFICATIONS

[sharp, toothed](#)