

## ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth

13.0cm (5.2")




**SKU:** 070-860-13

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV085
<b>KLS Martin</b>	15-770-13
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV085
<b>KLS Martin</b>	15-770-13
<b>AUTHORS</b>	<a href="#">ADSON</a>

## SPECIFICATIONS

[sharp, toothed](#)