

ADSON DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 2.0 +/-0.2


SKU: 020-070-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD223
KLS Martin	12-168-15
AUTHORS	ADSON , MICRO
Medicon	06.22.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD223
KLS Martin	12-168-15

AUTHORS

[ADSON](#), [MICRO](#)

Medicon

06.22.15

SPECIFICATIONS

[serrated](#), [straight](#)