

## ADSON DRESSING FORCEPS

straight; serrated

12.0cm (4.8"); width of jaw 1.3 +/-0.1


**SKU:** 020-070-12

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD222
<b>KLS Martin</b>	12-168-12
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	06.20.12
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD222
<b>KLS Martin</b>	12-168-12

## AUTHORS

[ADSON](#)

## Medicon

06.20.12

## SPECIFICATIONS

[serrated](#), [straight](#)