

ADSON HAEMOSTATIC FORCEPS

straight; serrated

18.5cm (7.4")

SKU: 040-414-18

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
Aesculap	BH200
KLS Martin	13-364-18
AUTHORS	ADSON
Medicon	16.40.19
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
Aesculap	BH200
KLS Martin	13-364-18

AUTHORS

[ADSON](#)

Medicon

16.40.19

SPECIFICATIONS

[serrated](#), [straight](#)