

## ADSON HAEMOSTATIC FORCEPS

straight; serrated

18.5cm (7.4")



**SKU:** 040-414-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH200
<b>KLS Martin</b>	13-364-18
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	16.40.19
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH200
<b>KLS Martin</b>	13-364-18

## AUTHORS

[ADSON](#)

## Medicon

16.40.19

## SPECIFICATIONS

[serrated](#), [straight](#)