

ADSON SELFRETAINING RETRACTOR

blunt; 3X4 teeth

28.0 cm; 11"




SKU: 070-683-26

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-737-26
AUTHORS	ADSON
Medicon	24.28.56
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-737-26
AUTHORS	ADSON
Medicon	24.28.56

SPECIFICATIONS

[blunt](#), [toothed](#)