

ADSON SELFRETAINING RETRACTOR

blunt; 6X6 teeth

29.0cm (11.6")




SKU: 070-683-29

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-739-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-739-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed