

ADSON SELFRETAINING RETRACTOR

long teeth left; sharp; 6X6 teeth

28.0cm (11")




SKU: 070-686-01

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV254
KLS Martin	15-790-02
AUTHORS	ADSON
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV254
KLS Martin	15-790-02
AUTHORS	ADSON

SPECIFICATIONS

[sharp, toothed](#)