

ADSON SELFRETAINING RETRACTOR

sharp; 3X4 teeth

28.0 cm; 11"




SKU: 070-682-26

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV252
KLS Martin	15-736-26
AUTHORS	ADSON
Medicon	24.28.55
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV252
KLS Martin	15-736-26

AUTHORS

[ADSON](#)

Medicon

24.28.55

SPECIFICATIONS

[sharp, toothed](#)