

## ADSON SUCTION TUBE

11 Fr.; angled

14.5cm (5.8") working length

**SKU:** 085-050-11

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">ADSON</a>         |

### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">ADSON</a>         |