

## ADSON SUCTION TUBE

15 Fr.; angled

14.5cm (5.8") working length

**SKU:** 085-050-15

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">ADSON</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">ADSON</a>