

ALLIGATOR FORCEPS

FL 600mm; 4 Fr.; double action

semi-rigid; with LuerLock; insert only



with Luer Lock

SKU: 583-460-12LLA

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments