

## ALLIGATOR FORCEPS; FL 600MM

5 Fr. semi-rigid; double action

with LuerLock; insert only



with Luer Lock

**SKU:** 583-460-16LLA

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscopy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscopy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>