

ANDERSON-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; curved

19.0cm (7.6")




SKU: 070-678-20

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV250
KLS Martin	15-732-20
AUTHORS	ADSON , ANDERSON
Medicon	24.28.02
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV250
KLS Martin	15-732-20

AUTHORS

[ADSON](#), [ANDERSON](#)

Medicon

24.28.02

SPECIFICATIONS

[sharp](#), [toothed](#)