

## ANDERSON-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; curved

19.0cm (7.6")




**SKU:** 070-678-20

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV250
<b>KLS Martin</b>	15-732-20
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">ANDERSON</a>
<b>Medicon</b>	24.28.02
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV250
<b>KLS Martin</b>	15-732-20

**AUTHORS**

[ADSON](#), [ANDERSON](#)

**Medicon**

24.28.02

**SPECIFICATIONS**

[sharp](#), [toothed](#)