

ANDREWS-PYNCHON SUCTION TUBE

removable suction tip

inner Ø 2.0mm; outer Ø 4.0mm; 24.0cm

SKU: 085-006-30

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	ANDREWS-PYNCHON
Medicon	09.85.08

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	ANDREWS-PYNCHON
Medicon	09.85.08