

## ANDREWS-PYNCHON SUCTION TUBE

removable suction tip

inner Ø 2.0mm; outer Ø 3.0mm; 24.0cm

**SKU:** 085-006-27

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">ANDREWS-PYNCHON</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">ANDREWS-PYNCHON</a>