

## MINI ENDOSCOPE HD+

Ø 1.9mm; 30°; 60 mm

autoclavable



**SKU:** 601-122-30HD

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Arthroscopy</a> , <a href="#">ENT</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Arthroscopy</a> , <a href="#">ENT</a>