

ARTHROSCOPE/ SINUSCOPE HD+

Ø 2.7mm; 70°; 187 mm

autoclavable


SKU: 601-203-70HD

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Arthroscopy , ENT |



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Arthroscopy , ENT |