

## ARTHROSCOPE/ SINUSCOPE HD+

Ø 4.0mm; 45°; 175mm

autoclavable

**SKU:** 601-204-40HD

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Arthroscopy</a> , <a href="#">ENT</a>



45°

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Arthroscopy</a> , <a href="#">ENT</a>