

ARTHROSCOPE/ SINUSCOPE HD+2

Ø 4.0mm; 30°; 175mm

autoclavable

SKU: 601-204-30HD2

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Arthroscopy](#), [ENT](#)