

ARTHROSCOPE/ SINUSCOPE HD+2

Ø 4.0mm; 45°; 175mm

autoclavable

SKU: 601-204-40HD2

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Arthroscopy , ENT |



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Arthroscopy , ENT |