

ARTHROSCOPY TELESCOPE HD

Ø 2.7mm; 70°; 110 mm

autoclavable

SKU: 601-303-70

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	Arthroscopy , ENT



70°

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	Arthroscopy , ENT