

## ARTHROSCOPY TELESCOPE HD

Ø 2.7mm; 0°; 110 mm

autoclavable

**SKU:** 601-303-00

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                                                   |
|--------------|---------------------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Telescopes</a>                        |
| <b>SCOPE</b> | <a href="#">Arthroscopy</a> , <a href="#">ENT</a> |



### ZUSÄTZLICHE INFORMATIONEN

|              |                                                   |
|--------------|---------------------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                              |
| <b>RANGE</b> | <a href="#">Telescopes</a>                        |
| <b>SCOPE</b> | <a href="#">Arthroscopy</a> , <a href="#">ENT</a> |