

## AUTOMATIC TROCAR -SLEEVE, Ø 13 MM,

WL120mm, stainless steel, maintenance free

stopcock, rough tube, distal bevelling



**SKU:** BLUO/04-00411-2-F

**Categories:** [Blue endo](#), [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)