

AUTOMATIC TROCAR -SLEEVE, Ø 5.5MM, WL75MM

short, Stainless Steel, maintenance-free

stopcock, rough tube, distal bevelling



SKU: BLUO/04-00415-2-SF

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)