

## BABY-ADSON HAEMOSTATIC FORCEPS

strongly curved; serrated

18.0cm (7.2")



**SKU:** 040-411-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ014
<b>KLS Martin</b>	13-421-18
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">BABY</a>
<b>Medicon</b>	17.05.19
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ014
<b>KLS Martin</b>	13-421-18

**AUTHORS**

[ADSON](#), [BABY](#)

**Medicon**

17.05.19

**SPECIFICATIONS**

[curved](#), [serrated](#)