

## BABY-DOYEN INTESTINAL CLAMP

straight; diagonally serrated

16.5cm (6.6")



**SKU:** 160-200-16

**Categories:** [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA120
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">DOYEN</a>
<b>SPECIFICATIONS</b>	<a href="#">diagonally serrated</a> , <a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA120
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">DOYEN</a>
<b>SPECIFICATIONS</b>	<a href="#">diagonally serrated</a> , <a href="#">straight</a>