

BABY-MIKULICZ PERITONEAL CLAMP

curved; serrated; 1X2 teeth

14.0cm (5.6")




SKU: 040-701-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Peritoneal Clamp Forceps |
| Aesculap | BJ300 |
| KLS Martin | 13-721-14 |
| AUTHORS | BABY , MIKULICZ |
| Medicon | 16.85.14 |
| SPECIFICATIONS | curved , serrated , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Peritoneal Clamp Forceps |
| Aesculap | BJ300 |
| KLS Martin | 13-721-14 |

AUTHORS

[BABY](#), [MIKULICZ](#)

Medicon

16.85.14

SPECIFICATIONS

[curved](#), [serrated](#), [toothed](#)