

BABY-MIXTER DISSECTING FORCEPES

strongly curved; serrated

14.0cm (5.6")




SKU: 040-475-02

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ008
KLS Martin	13-303-02
AUTHORS	BABY , MIXTER
Medicon	17.06.14
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ008
KLS Martin	13-303-02

AUTHORS

[BABY](#), [MIXTER](#)

Medicon

17.06.14

SPECIFICATIONS

[curved](#), [serrated](#)