

## BABY-MIXTER DISSECTING FORCEPS

slightly curved; serrated

14.0cm (5.6")




**SKU:** 040-475-01

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ011
<b>KLS Martin</b>	13-471-14
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">MIXTER</a>
<b>Medicon</b>	17.06.14
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ011
<b>KLS Martin</b>	13-471-14

## AUTHORS

[BABY](#), [MIXTER](#)

## Medicon

17.06.14

## SPECIFICATIONS

[curved](#), [serrated](#)